medicura®

ISTRUZIONI PER L'USO
INSTRUCTIONS FOR USE
FOLLETO DE INSTRUCCIONES

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Andon Health Co.,LTD
No.3 Jinping Street, Ya An Road, Nankai District, Tianjin 300190, China

Lotus Global Co., LTD - 1 Four Seasons Terrace, West Drayton
Middlesex, London, UB7 9GG, United Kingdom
Tel.:+0044-20-75868010 - Fax: +0044-20-79006187

WARNINGS

ATTENTION : read this booklet carefully it contains important instructions for a safe installation, use and maintenance.

Important instructions to be kept for future reference.

Keep handy the guarantee, receipt of payment and these instructions for future use and reference.

After removing the packaging, check the appliance for any signs of damage.

In case of visible damage, do not use it and contact a qualified technician.

Keep all parts of the packaging out of the reach of children.

The appliance is not a toy: it is an electrical appliance and must be treated with the necessary caution.

The use of any electrical appliance requires to respect basic rules, namely:

- Never touch the appliance with wet or humid hands.

- his appliance must not be used by individuals.

(including children) with diminished physical, sensory or mental abilities or who lack the experience or knowledge on how to operate it unless they are supervised or have been instructed on the use of the appliance by someone who is responsible for their safety. Children must be supervised to make sure that they do not play with the appliance.

- Do not leave the appliance exposed to weathering (rain, sun).

- Do not subject it to rough handling.

It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer. Please do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.

Please do not share the cuff with other infective person to avoid cross-infection.

If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of will be displayed.

Under this condition, the Electronic Sphygmomanometers can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment. In case of a breakdown or malfunctioning of the appliance, do not attempt to repair it, but take it to a service centre approved by the manufacturer. For safety reasons, the appliance should never be opened.


RESPONSIBILITIES

The manufacturer, assembler, installer or importer will be responsible for safety, performance and reliability only if:

- The assembly, calibration, repair and modification have been carried out by duly authorised persons;
- The instructions have been followed;
- The manufacturer takes no responsibility for damage caused by incorrect, improper or unreasonable use or the incorrect use of the accessories.

MEASUREMENT WARNINGS


Avoid eating, smoking or performing intense activities at least 30 minutes before taking your blood pressure. Rest for about 5/10 minutes before taking your blood pressure. You need to be seated when taking your blood pressure. Make sure that the arm cuff is at the same height as the heart and properly wrapped around the arm. If desired, you can lie down, always making sure that your arms are stretched out and the palm is open and turned upwards. Make sure that the cuff is at the same level as the heart.

Do not speak or move during measurement. Take your blood pressure always on the same arm, failure to do so may account for significant differences. The time between two measurements must be 3 minutes to allow for circulation recovery. The cuff fits normal adult arm size. Measurement is not possible for people who suffer cardiac arrhythmia. Avoid electromagnetic interferences caused by mobile phones, electric razors, microwave ovens, etc. Wrap the cuff directly on the skin for accurate blood pressure and pulse measurements.

TECHNICAL INFORMATIONS.

- 1 - Display
- 2 - Pressure level indicator
- 3 - On/Off key
- 4 - Memory key
- 5 - Mode / date/time setting key
- 6 - + / - keys
- 7 - Voice function key

Model: m245/KD-5904.

Measurement method: oscillometric, automatic inflation and measurement. Storage capacity: 60 sets of measurement.

Arm cuff circumference: 22 ~ 30 cm

Appliance dimensions: 145x104x55 mm.

Weight: 310 g (without batteries).

Battery life : approx 360 times.

Accuracy of blood pressure measurement: ± 3 mmHg.

Blood pressure measurement tolerance: 0-300 mmHg.

Heart beat measurement tolerance: 40-180/min

Alimentation: batteries 4x1,5V  type AA or DC 6V 600mA (adaptor not included).

Transport and storage conditions.

Humidity tolerance: $\leq 90\%$ RH.

Room temperature tolerance: $-20 \sim 55^{\circ}\text{C}$.

Functioning conditions.

Humidity tolerance: $\leq 90\%$ RH.

Room temperature tolerance: $5 \sim 40^{\circ}\text{C}$.

The list of all components belonging to the pressure measuring system, including accessories:

Pump, Valve, LCD, Cuff, Sensor.

SETTING THE CLOCK

With the appliance off, press the MODE key; the year will start flashing and can be set using the +/- keys.

Press the MODE key again to select the day and then the time. To adjust the settings, press the +/- keys.

CAUTION: The date and time must be set again after replacing the batteries!

VOICE FUNCTION SETTING

With the appliance off, press the voice function key and hold it pressed. Without releasing the key, view the available languages scrolling on the display (L1, L2, etc.).

When the desired language appears, release the key.

Selecting "L0" will switch off the voice function.

Then press the voice function key again to set the voice volume. The volume will increase each time the key is pressed.

INSTRUCTION FOR USE

Insert the rubber tube into the apparatus, making sure that no air escapes. Place the arm band over the left arm, positioning it approximately 1-2 cm above the elbow. Keep the palm of the hand facing upward and move the rubber tube so that it is parallel to the arm and positioned on the artery. Make sure that the arm cuff is not too tight (you should be able to slip a finger between the arm cuff and your arm).

Measurement

Press the "START" button to activate the appliance. All the symbols on the display will come on, after which the last measurement performed will appear. After a few seconds the arm cuff will inflate automatically.

The arm cuff will be inflated to the pressure of about 190mmHg and then slowly deflate while measuring the pressure of the subject.


At this point the arm cuff deflates and the heart symbol flashes on the display, in synchrony with the heart beat, if the pulses can be detected.

Relax and do not move during the measurement.

At the end of the measurement, the systolic pressure (maximum), diastolic pressure (minimum), pulse rate and pressure level indicator appear on the display.

If the voice function has been activated, the appliance will say the values recorded and the classification level in the preselected language.

If the appliance records an irregular pulse value, the "HEART" symbol appears on the display.

The appliance automatically stores the result for up to a maximum of 60 measurements. When this limit has been reached, the new results progressively replace the previous ones. The stored data can be displayed complete with date and time.

During measurement, you can press the “START” button to turn off the monitor manually.

Repeating the measurement.

Switch the appliance off, then simply press the “START” push button to start another measurement.

Wait at least three minutes before the next measurement.

Calling up stored results.

With the appliance off, press the MEM key to start the memory data call-up mode.



Each time the MEM key is pressed, the recorded measurements are displayed.

If the voice function has been activated, the appliance will say the memory position number and the values stored.

If no measurement has been recorded, “00” will appear on the display as the memory position number and “0” will appear as the recorded value.

To cancel the saved data, press the MEM key and hold it pressed for about 3 seconds until the appliance emits three beeps.

The display will appear first as shown in fig. 4 and then as shown in fig. 5.

Battery installation or replacement

When the battery symbol appears on the display, this means that the batteries must be replaced. Proceed following the instructions below.

Make certain that the appliance is turned off.

Remove the cover and insert 4x1,5V AA alkaline batteries, according to the polarity symbols. Replace the cover correctly.

If the appliance will not be used for an extended period of time, it is recommended that you remove the batteries.

MAINTENANCE.

This appliance does not require special maintenance. Clean it regularly.

Use a soft, damp cloth to clean the appliance and the cuff.

Air dry the cuff and do not iron it.

Do not use abrasive products or chemical solvents.

DO NOT IMMERSER THE APPLIANCE IN WATER.

No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairably can be supplied.

It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.

Avoid high temperature and solarization

It is recommended the cuff should be disinfected 2 times every week if needed (For example, in hospital or in clinique). Wipe the inner side (the side contacts skin) of the cuff by a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing. 2. Clean the cuff after the usage of every 200 times is recommended.

The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years

The cuff integrity is maintained after 1,000 open-close cycles of the closure.

Classification of the measurements

Classification of the values	Systolic Pressure	Diastolic Pressure	Indicator colour
Excellent	<120	<80	green
Normal	120-129	80-84	green
Normal-High	130-139	85-89	green
Hypertension Level 1	140-159	90-99	yellow
Hypertension Level 2	160-179	100-109	orange
Hypertension Level 3	> 180	> 110	red

TROUBLESHOOTING

Type of error	Cause	Checks
Battery symbol appears on the	down batteries	Do not move, and repeat the measurement
Er0	The user has moved during the measurement	Do not move, and repeat the measurement
Er1	Error in the measurement of the systolic pressure	Repeat the measurement
Er2	Error in the measurement of the diastolic pressure	Repeat the measurement
Er3	Inflation blocked because the wristband is too tight	Reposition the wristband correctly and check the connection
Er4	Air leak or the wristband is too loose	Reposition the wristband correctly and check the connection
Er5	Inflation above 300 mmHg	Press the START key
Er6	Inflation above 15 mmHg after 3 minutes	Press the START key

Er7	Error in accessing the memory system	Contact the retailer
The entire display is flashing	Error in the control of the appliance parameters	Contact the retailer
Er 9	Error in the self-diagnosis of the MCU	Contact the retailer
ErA	Capacity error	Contact the retailer
ErB	Error when saving in the memory system	Contact the retailer
The appliance does not respond when the keys are pressed	Incorrect operation	Remove the batteries for 5 minutes and then reinstall them

Symbol for "THE OPERATION GUIDE MUST BE READ"
(The sign background colour: blue.The sign graphical symbol: white)

Symbol for "WARNING"

Symbol for "TYPE BF APPLIED PARTS"
(The cuff is type BF applied part)

Symbol for "MANUFACTURER"

CE 0197

Symbol for "COMPILES WITH MDD93/42/EEC REQUIREMENTS"

Symbol for "DATE OF MANUFACTURE"

Symbol for "EUROPEAN REPRESENTATION"

SN

Symbol for "SERIAL NUMBER"

Symbol for "KEEP DRY"

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

m245 is intended for use in the electromagnetic environment specified below.
The customer or the user of the m245 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The m245 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The m245 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity			
The m245 is intended for use in the electromagnetic environment specified below. The customer or the user of the m245 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity

The m245 is intended for use in the electromagnetic environment specified below.
 The customer or the user of the m245 should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the m245, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3\sqrt{P} \text{ 80 MHz to 2,5 GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the

			recommended separation distance in meters (m).
--	--	--	--

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the m245 is used exceeds the applicable RF compliance level above, the m245 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the m245.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

**Recommended separation distances between
portable and mobile RF communications equipment and the m245**

The m245 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the m245 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the m245 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

EN 60601-1:2006/AC:2010 (Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance), EN 60601-1-2:2007 /AC:2010 (Medical electrical equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility -Requirements and tests), EN 1060-1: 1995 + A1:2002 + A2: 2009 (Noninvasive sphygmomanometers - Part 1: General requirements), EN 1060-3: 1997 + A1:2005 + A2: 2009 (Noninvasive sphygmomanometers - Part 3: Supplementary requirements for electromechanical blood pressure measuring systems), ANSI/AAMI SP-10:2002+A1:2003+A2:2006.

Andon Health Co.,LTD
No.3 Jinping Street, Ya An Road, Nankai District, Tianjin 300190, China

Lotus Global Co., LTD - 1 Four Seasons Terrace, West Drayton
Middlesex, London, UB7 9GG, United Kingdom
Tel.:+0044-20-75868010 - Fax: +0044-20-79006187

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Terms and Conditions

The guarantee is valid for 24 months.

This guarantee applies only if it has been duly filled in and is submitted with the receipt showing the date of purchase.

The guarantee covers the replacement or repair of parts making up the appliance which were faulty at source due to manufacturing faults.

After the guarantee has expired, the appliance will be repaired against payment.

The manufacturer declines any responsibility for damage to persons, animals or property due to misuse of the appliance and failure to observe the directions contained in the instructions.

Limits

All rights under this guarantee and any responsibility on our part will be voided if the appliance has been:

- mishandled by unauthorized persons
- improperly used, stored or transported.

The guarantee does not cover damage to the outward appearance or any other that does not prevent regular operation.

If any faults should be found despite the care taken in selecting the materials and in creating the product, or if any information or advice are required, please contact your local dealer.

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC.

At the end of its working life, the product must not be disposed of as urban waste.

It must be taken to a special local authority differentiated waste collection centre or to a dealer providing this service.

Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant savings in energy and resources. As a reminder of the need to dispose of household appliances separately, the product is marked with a crossed-out wheeled dustbin.

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Prodotto importato da POLY POOL S.p.A.